



*Quality you trust. Experience you know.*

## *Xebeco Ground Beef Fat Analyzer*



### *Features:*

- The industry standard for measuring the fat content of ground beef.*
- Fast, simple, and affordable*
- Single lever manual operation takes only seconds*
- Three minute cycle time*
- Standard voltage:  
115V/60Hz/1Ph*
- Comes complete with carrying case and glass vials*
- The XeFA is a self calibrating machine as each time you measure the fat content you are setting the machine to read zero. No external calibration is needed.*

### *Models:*

*XeFA*

**MADE IN THE USA**

### *Standard Accessories:*

- Sample measuring device*
- Carrying case*
- 6 glass vials*

### *Optional Accessories:*

- 1400035 Carton of 48 glass vials*


### *Warranty:*

- 90 day, carry-in parts and labor*

#### *Warranty Details*

*The Fat Analyzer units carry a 90 day, carry-in parts and labor warranty against any flaws in materials and workmanship. This period begins on the date of the purchase and, provided the unit is used in accordance with our instructions, is in full effect for 90 days from that date. Any work performed under this warranty must be performed between the hours of 8am and 5pm Monday through Friday. Xebeco will not pay for overtime charges for work performed other than during normal business hours. Damage incurred in transit or from installation error, accident, alteration or misuse is not covered. Xebeco shall not be liable for any consequential compensatory, incidental, or special damages. Contact the Xebeco service department for any warranty claims.*

Quality you trust. Experience you know.



A: Width  
5.5" (140mm)

B: Depth  
5.4" (137mm)

C: Height  
5.6" (142mm)

**Electrical:** 115V/60Hz/1PH

**Finish:** Black wrinkle enamel.

**Calibration:** The XeFA is a self calibrating machine as each time you measure the fat content you are setting the machine to read zero. No external calibration is needed.

**Standard Equipment:** Sample measuring device, carrying case, and 6 glass vials.

**Optional Accessories:** Case of 48 additional glass vials.

**Accuracy:** The XeFA was tested against the AOAC method of determining fat content as well as the Modified Babcock method. The test was done with three samples of regular, lean, and extra lean meat being tested and averaged; the results are within 1%.

	AOAC	Modified Babcock	Xebeco
Regular	23.84	24.7	24.1
Lean	21.58	21.8	22.1
Extra Lean	15.21	16.1	16.2

### Shipping Information:

Freight Class: 85  
Net Weight: 8 lbs (3.6Kg)  
Shipping Weight: 10 lbs (4.5Kg)

### Shipping Box Dimensions:

Height: 12" (305mm)  
Width: 12" (305mm)  
Depth: 12" (305mm)  
Cubic Measure: 1ft<sup>3</sup> (0.03M)<sup>3</sup>

Cartons are affixed to a pallet. The pallets are heat-treated and NELMA approved for international shipping. The weight and dimensions of the pallet are included in the weight and dimensions above.

MADE IN THE USA